

Configuração Completa - LERP Intranet

18 Categorias de Reagentes + 10 Localizações



CATEGORIAS DE REAGENTES (18)

#	Categoria	Letra	Código Interno	Local de Armazenamento Padrão
1	Flammable	U	LERP-U####	CFF - Flammable (all categories)
2	Acid	A	LERP-A####	CAB 1 — SHF 1 - Acid
3	Corrosive	U	LERP-U####	CAB 1 — SHF 2 - Corrosive (non-acid, non-base)
4	Base	B	LERP-B####	CAB 2 — SHF 1 - Base
5	Inorganic Salt	I	LERP-I####	CAB 2 — SHF 2 - Inorganic Salt
6	Polymer	P	LERP-P####	CAB 3 — SHF 1 - Polymer/Thiol/Catalyst/Cross-linker/Photoinitiator/Solvent (inert)
7	Thiol	T	LERP-T####	CAB 3 — SHF 1 - Polymer/Thiol/Catalyst/Cross-linker/Photoinitiator/Solvent (inert)
8	Catalyst	C	LERP-C####	CAB 3 — SHF 1 - Polymer/Thiol/Catalyst/Cross-linker/Photoinitiator/Solvent (inert)
9	Crosslinker	X	LERP-X####	CAB 3 — SHF 1 - Polymer/Thiol/Catalyst/Cross-linker/Photoinitiator/Solvent (inert)

#	Categoria	Letra	Código Interno	Local de Armazenamento Padrão
10	Photoinitiator	F	LERP-F####	CAB 3 — SHF 1 - Polymer/Thiol/ Catalyst/Cross-linker/Photoinitiator/Solvent (inert)
11	Solvent (inert)	S	LERP-S####	CAB 3 — SHF 1 - Polymer/Thiol/ Catalyst/Cross-linker/Photoinitiator/Solvent (inert)
12	Oxidizer / Reducer	O	LERP-O####	CAB 3 — SHF 2 - Oxidizer/Nanomaterial/ Analytical/ Monomer
13	Nanomaterial	N	LERP-N####	CAB 3 — SHF 2 - Oxidizer/Nanomaterial/ Analytical/ Monomer
14	Analytical	L	LERP-L####	CAB 3 — SHF 2 - Oxidizer/Nanomaterial/ Analytical/ Monomer
15	Monomer	M	LERP-M####	CAB 3 — SHF 2 - Oxidizer/Nanomaterial/ Analytical/ Monomer
16	Controlled	K	LERP-K####	CAB 4 — SHF 1/2 - Controlled
17	Microbiology	G	LERP-G####	CAB 5 — SHF 1/2 - Microbiology

#	Categoria	Letra	Código Interno	Local de Armazenamento Padrão
18	Refrigerated	U	LERP-U####	Refrigerator/ Freezer - Refrigerated

Nota: #### = 4 dígitos aleatórios gerados automaticamente (ex: LERP-A1234)

LOCALIZAÇÕES DE ARMAZENAMENTO (10)

#	Nome	Descrição	Tipo
1	CFF	Flammable Storage Cabinet	Fume Hood
2	CAB 1 — SHF 1	Cabinet 1, Shelf 1 — Acids	Cabinet
3	CAB 1 — SHF 2	Cabinet 1, Shelf 2 — Corrosive	Cabinet
4	CAB 2 — SHF 1	Cabinet 2, Shelf 1 — Bases	Cabinet
5	CAB 2 — SHF 2	Cabinet 2, Shelf 2 — Inorganic Salts	Cabinet
6	CAB 3 — SHF 1	Cabinet 3, Shelf 1 — Polymers/Catalysts/ etc.	Cabinet
7	CAB 3 — SHF 2	Cabinet 3, Shelf 2 — Oxidizers/Monomers/ etc.	Cabinet
8	CAB 4 — SHF 1/2	Cabinet 4 — Controlled Substances	Cabinet
9	CAB 5 — SHF 1/2	Cabinet 5 — Microbiology	Cabinet
10	Refrigerator/Freezer	Laboratory Refrigerator / Freezer	Refrigerator

MAPEAMENTO POR ARMÁRIO

CFF (Capela de Exaustão)

- Flammable

CAB 1 — Armário 1

- SHF 1 (Prateleira 1): Acid
- SHF 2 (Prateleira 2): Corrosive

CAB 2 — Armário 2

- SHF 1 (Prateleira 1): Base
- SHF 2 (Prateleira 2): Inorganic Salt

CAB 3 — Armário 3

- SHF 1 (Prateleira 1): Polymer, Thiol, Catalyst, Crosslinker, Photoinitiator, Solvent (inert)
- SHF 2 (Prateleira 2): Oxidizer / Reducer, Nanomaterial, Analytical, Monomer

CAB 4 — Armário 4

- SHF 1/2: Controlled

CAB 5 — Armário 5

- SHF 1/2: Microbiology

Refrigerator/Freezer

- Refrigerated
-

Como Gerenciar

Acesse: **Admin** → **Categories & Storage Settings** (</admin/settings>)

- **Adicionar/Editar/Remover** categorias e localizações dinamicamente
 - O campo “**Letter**” define a letra do código interno (LERP-**{letter}**{4digits})
 - O campo “**Default Storage Location**” define o preenchimento automático no formulário
-

Re-executar Seed

Para restaurar todos os dados originais:

```
npx tsx prisma/seed-complete-data.ts
```

Data: Maio 2026

Sistema: LERP Intranet v2

Laboratório: LERP — FEQ/UNICAMP